

# P.T.Rail™ Floor Mast

Toilet Support Rail Floor Mast

## Technical Data Sheet

### Product Specification

The grab rail floor support mast shall be specifically designed for the HealthCraft PT Rail. The product shall be compatible with either a fold up grab rail or a fixed grab rail. The product shall allow the top of the grab rail to be mounted at a height ranging from 31"/78.7 cm to 36"/91.4 cm. The total height of the product shall be 39"/99.1 cm. The base plate of the product shall be 0.25"/0.64cm thick and have 5 fastener holes. The product shall have a weight capacity of 400 lbs/181 kg. The grab rail floor support mast shall be HealthCraft Products model PT-WRM-\_\_\_\_ (see table below for appropriate product code).



### Key Features

- Compatible with Both Fixed and Folding Rails
- Flexible Mounting Height
- Solid Construction
- 400bs / 181kg Weight Capacity

### Codes and Standards

- ADA Compliant Structural Strength (ADAAG 4.26.3)

### Specified Model


Model	Finish	Material
PT-WRM-PCW	Powder Coat White	Steel
PT-WRM-SS	Stainless Steel	Stainless Steel

# P.T.Rail™ Floor Mast

Toilet Support Rail Floor Mast

## Technical Data Sheet

### Product Dimensions



Model	Height (H)	Width (W)	Depth (D)	Weight
PT-WRM	6" / 15.2cm	39" / 99.1cm	7" / 17.8cm	22lbs / 10.0kg

