


Technical Data Sheet

Product Specification

The support rail shall comprise 2 rails: a fixed bed rail and a pivot rail. The fixed rail shall be 36" / 91cm in length, installed between the mattress and box spring and secured in place with a safety strap. The fixed horizontal handle shall have an overall 1.5" / 3.8cm rail diameter, extending 19.75" / 50cm across. A storage compartment shall span the opening underneath the handle. The fixed bed rail shall be height adjustable between 18-22" / 46-56cm. The pivot rail shall have an overall 1.5" / 3.8cm rail diameter, extending 16.25" / 41cm across. It shall contact the ground and pivot freely 180°. The pivot rail shall be height adjustable between 29.25-38.25" / 74-97cm. The product shall have a user weight capacity of 300 lbs/136 kg, and is intended only for moderate vertical load bearing. The support rail shall be HealthCraft Products model SRS-S.



Key Features

- 180° Pivoting Support Rail
- Comfortable EVA Foam Grip
- Storage Compartment
- Height Adjustable



Codes and Standards


Complies with ASTM F3186-17 - Standard Specification for Adult Portable Bed Rails and Related Products.

Specified Model

Model	Weight Capacity	Finish	Material
SRS-S	300 lbs / 136 kg	White Gloss	Steel

Technical Data Sheet

Product Dimensions



Shipping Dimensions

Part	Height (H)	Width (W)	Depth (D)	Weight
SRS-S	3.3" / 8.5cm	17" / 43cm	28" / 71cm	12.40lbs / 5.6kg

